



a passion for innovation®

www.ldrmedical.com

France
Hôtel de Bureaux 1,
4 rue Gustave Eiffel,
10430 Rosières Près Troyes
Mailing Address:
10902 Troyes Cedex 9, France
Technopôle de l'Aube BP 2
+33 (0)3 25 82 32 63

China
Beijing Global Trade Center # 36
North Third Ring Road East,
Unit 06, Level 19, Building A,
Dongcheng District,
Beijing, China, 100013
+86 10 58256655

Brazil
Av. Pereira Barreto, 1395
Torre sul - CJ 193 - Bairro Paraiso
Santo André - São Paulo
CEP: 09190-610
Brazil
+55 11 43327755

United States
13785 Research Boulevard
Suite 200
Austin, TX 78750
USA
512.344.3333



Document intended for the exclusive use of healthcare professionals.
C-Plate® - Non-sterile anterior cervical plate system - is a class IIb CE marked medical device made by the LDR Medical S.A.S. Company and for which the conformity assessment was carried out by the notified body G-Med N°0459. C-Plate® is an anterior cervical plate system intended to provide immobilization and stabilization of spinal segments during the development of a cervical spinal fusion. Before any surgical procedure, read carefully the instructions and the surgical technique.

LDR, LDR Spine, LDR Medical, BF+, BF+(ph), EasySpine, C-Plate, SpineTune, MC+, Mobi+, Mobic, Mob-L, Mobicis, ROI, ROI-A, ROI-T, ROI-C, Avenue L and verteBRIDGE are trademarks or registered trademarks of LDR Holding Corporation or its affiliates in France, the United States or other countries.



a passion for innovation®

C-Plate®

ANTERIOR CERVICAL PLATE




C-Plate®

ANTERIOR CERVICAL PLATE



A secured assembly




A complete range

A full range of **semi-fixed** and **variable screws** and **plates** cover multiple indications while filling the surgeons needs.




A zero profile design

Only **1,6 mm thick**, C Plate is one of the lowest profile plating systems on the market.

Its narrow width, its atraumatic integrated locking mechanism and its anatomical design provides the patient added comfort and limits risks of post-surgery complications.




A streamlined procedure



From 10 to 18 mm length



Screw insertion angulation



- All C-Plates are delivered **pre-formed** (sagittal and axial plans) to precisely fit the cervical spine while reducing operative time.
- Tapping is optional thanks to the use of **self-tapping** screws.
- Screw prehension from the implant rack is "no touch" due to the self-retaining print of the screw heads.

